

Test Items

B-CUP 20



Leading brand in Vitro Diagnostics

Cat. No.	Product Description	Specimen	Format	Kit Size
Reader				
AOABCUP20	Cup Reader	/	Analyzer	1 unit
Rapid Tests for Reader				
AOA-R127-A1	Multi-Drug 2 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R137-A1	Multi-Drug 3 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R147-A1	Multi-Drug 4 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R157-A1	Multi-Drug 5 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R167-A1	Multi-Drug 6 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R177-A1	Multi-Drug 7 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R187-A1	Multi-Drug 8 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R197-A1	Multi-Drug 9 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1107-A1	Multi-Drug 10 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1117-A1	Multi-Drug 11 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1127-A1	Multi-Drug 12 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1137-A1	Multi-Drug 13 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1147-A1	Multi-Drug 14 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1157-A1	Multi-Drug 15 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1167-A1	Multi-Drug 16 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1177-A1	Multi-Drug 17 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AOA-R1187-A1	Multi-Drug 18 Drugs Rapid Test 1-Step Cup	Urine	Cup	25T
AUA-R127-A1	Multi-Drug 2 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R137-A1	Multi-Drug 3 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R147-A1	Multi-Drug 4 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R157-A1	Multi-Drug 5 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R167-A1	Multi-Drug 6 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R177-A1	Multi-Drug 7 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R187-A1	Multi-Drug 8 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R197-A1	Multi-Drug 9 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1107-A1	Multi-Drug 10 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1117-A1	Multi-Drug 11 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1127-A1	Multi-Drug 12 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1137-A1	Multi-Drug 13 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1147-A1	Multi-Drug 14 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1157-A1	Multi-Drug 15 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1167-A1	Multi-Drug 16 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T
AUA-R1177-A1	Multi-Drug 17 Drugs Rapid Test 1-Step Cup with Adulteration	Urine	Cup	25T

B-CUP 20

• Cup Reader •

PREMIUM POINT-OF-CARE SOLUTION



DrugLine™



Biozek B.V.



Laan van de Ram 49

7324BW, Apeldoorn The Netherlands



info@biozek.com



www.biozek.com

Advantages

User-Friendly

7" touch screen
Fast result at 5 min
Automatic print
Easy operation

Portable Design

Small footprint
Lightweight


LIS Compatibility

- Anti-narcotics departments
- Public security bureaus
- Drug rehabilitation centers
- Physical examination centers
- National health inspection agencies
- Clinical Research Organizations (CROs)
- Employment screening
- Tourism companies

Technical Specifications

Principle	Photoelectric Conversion
Test Format	Cup
Touch Screen	7" LCD Touch Screen
Measurement	Qualitative
Read Time	1 - 2.5 Minutes
Power Source	DC 12V, 2A
Dimensions (L*W*H)	221* 215 *165 mm
Weight	< 1.5 kg
Connecting Ports	LAN Port, 4 USB Ports
Operating Conditions	15 - 30 °C
Printer	Built-in thermal printer

B-CUP 20 cup reader



Test Items

Multi-Drug Rapid Test 1-Step Cup With/Without Adulteration (Specimen:Urine)		
AB-PINACA (ABP)	Ecstasy (MDMA)	Phencyclidine (PCP)
Synthetic Marijuana (K2)	Tenamfetamine (MDA)	Pregabalin (PGB)
UR-144	Methylenedioxypropylvalerone (MDPV)	Propoxyphene (PPX)
Acetaminophen (ACE)	Methamphetamine (MET)	Quetiapine (QTP)
7-Aminoclonazepam (7-ACL)	Mephedrone HCl (MEP)	Risperidone (RPD)
Alcohol (ALC)	Mescaline (MES)	Scopolamine (SCOP)
Alprazolam (ALP)	Mirtazapine (MIT)	Tapentadol (TAP)
Amphetamine (AMP)	Morphine (MOP)	Tricyclic Antidepressants (TCA)
α-Pyrrolidinovalerophenone (α-PVP)	Methylphenidate(MPD)	Marijuana (THC)
Barbiturate (BAR)	Meperidine (MPRD)	Tilidine (TLD)
Buprenorphine (BUP)	Methaqualone (MQL)	Tramadol (TML)
Benzodiazepines (BZO)	Methadone (MTD)	Tropicamide (TRO)
Cathine (CAT)	N,N-Dimethyltryptamine (NND)	Trazodone (TZD)
Caffeine (CAF)	Opiates (OPI)	Zaleplon (ZAL)
Carisoprodol (CAR)	Oxycodone (OXY)	Zolpidem (ZOL)
Citalopram(CIT)	Olanzapine (OZP)	Zopiclone (ZOP)
Cannabinol (CNB)	Papaverine (PAP)	Urine Adulteration (OX/NIT/SG/GLU/pH/CRE)
Carfentanyl (CFYL)		
Clonazepam (CLO)		
Cocaine (COC)		
Cotinine (COT)		
Diazepam (DIA)		
Ethylenediamine-dimethylphosphinic acid (EDDP)		
Ethyl Glucuronide (ETG)		
Fluoketamine (FKET)		
Fluoxetine (FLX)		
Fentanyl (FYL)		
Gabapentin (GAB)		
γ-Hydroxybutyric acid (GHB)		
Ketamine (KET)		
Kratom (KRA)		
Lysergic Acid Diethylamide (LSD)		
6-Monoacetylmorphine (6-MAM)		
Methcathinone (MCAT)		



Convenient Operation

